




ORTH ORIGINAL

EVOCABLE™

Cerclage System



Periprosthetic fractures



Prophylactic cerclage



Greater trochanter reattachment

Multifilament design

High tensile strength

Superior fatigue properties




www.orthoriginal.bzh



Standard cable



Beaded cable



Sleeve



Available in two settings:

Standard Evocable set (sterile)

Part No : 0073-00-01

- Standard cable x1
- Sleeve x1

Standard cable

- Ø2,0mm cable, length 600mm
- In surgical stainless steel

Beaded cable

- Ø2,0mm cable, length 600mm, Ø3,0mm bead
- In surgical stainless steel

Sleeve

- Trapezoid shape
- In surgical stainless steel

Beaded Evocable set (sterile)

Part No : 0073-00-02

- Beaded cable x1
- Sleeve x1

For more information :

Distributed by

AKED!

44C rue de Bray
35510 Cesson-Sévigné,
FRANCE

Tél : +33 (0)2 23 35 53 40
Fax : +33 (0)2 23 35 16 35
contact@aked.bzh



Full instrument set

CE
0459