


Broken screw extraction set EMER Set



A simple and compact solution for
extraction of broken or damaged
screws
from $\varnothing 1.5\text{mm}$ to $\varnothing 7.0\text{mm}$



Conical extraction tips (mm)
 $\varnothing 1,5$ à $\varnothing 2,0$ - $\varnothing 2,7$ à $\varnothing 4,0$ - $\varnothing 4,5$ à $\varnothing 7,0$



AO coupling T-handle



Screw shaft removal tips (mm)

- $\varnothing 1,5$
- $\varnothing 2,0$
- $\varnothing 2,7$
- $\varnothing 3,5 - \varnothing 4,0$
- $\varnothing 4,5$
- $\varnothing 5,0 - \varnothing 7,0$



Hollow reamers with optional guide rod (mm)

- $\varnothing 1,5$
- $\varnothing 2,0$
- $\varnothing 2,7$
- $\varnothing 3,5 - \varnothing 4,0$
- $\varnothing 4,5$
- $\varnothing 5,0 - \varnothing 7,0$



Hook - Length 155mm



Spare Hollow reamers (mm)

- $\varnothing 1,5$
- $\varnothing 2,0$
- $\varnothing 2,7$
- $\varnothing 3,5 - \varnothing 4,0$
- $\varnothing 4,5$
- $\varnothing 5,0 - \varnothing 7,0$



Chisels
Width 4, 6, 10mm



Screw holding forceps - Length 205mm