



ORTHOORIGINAL

Graduated chisels

For tibial osteotomy




Open and closed-wedge
tibial osteotomy



Incremental width from
10mm to 25mm



Graduated, with a strike
plate



Tibial osteotomy chisels


3mm thickness

Strike plate

Sharp tip

10mm wide

REF 0092 01 10




15mm wide

REF 0092 01 15



20mm wide

REF 0092 01 20



25mm wide

REF 0092 01 25



240mm long


Graduated from 0 to 110mm

Can be stacked together

In stainless steel

C marked medical device

More sizes on request



For more information :



Tel. : +33 (0)2.23.35.53.40 contact@orthoriginal.bzh

Fax : +33 (0)2.23.35.16.35 www.orthoriginal.bzh

ORTHORIGINAL, 44C Rue de Bray
35510 Cesson-Sévigné, France