

Acetabular cup extraction system

To quickly remove an acetabular cup with minimal loss of bone



Hip revision surgery



For uncemented stems




Bone stock preservation

Optimized
number of
component

Quick and
easy to
assemble


Pre, per and
post
operative
time gain



Acetabular cup extraction system

Spherical heads


Half-spherical heads



Starter and finish blades

Universal handle

Cup caliper



Perforated box with lid



For a cup with ...



... or without insert.

Spherical heads

- For use in an acetabular cup where the insert is not removed (ceramic insert for example)
- Sizes : Ø22, Ø28, Ø32, Ø36 mm

Half-spherical heads

- For use in an acetabular cup where the insert has been removed.
- Sizes : from Ø33 to Ø49 (2mm increment)

Removable blades

- Two blades (starter and finish) per size
- Sizes : from Ø44 to Ø60 (2mm increment)

Handle

- Universal handle where every blades and heads fit.

Cup caliper

- To measure the diameter of head or the acetabular cup to remove

Contact us for
the full list of
part numbers

For more information:



Tél. : +33 (0)2.23.35.53.40 contact@orthoriginal.com
Fax : +33 (0)2.23.35.16.35 www.orthoriginal.com

Orthoriginal, 44C Rue de Bray
35510 Cesson-Sévigné, France

