




ORTH ORIGINAL

# Femoral stem extractor

Extraction of femoral stems made easy




Usable with all stems with a trunnion



Supplied with a pivoting slap hammer



Robust & reliable device



1. The extractor is compatible with all stems with a modular head, after removal of the head.



2. The trunnion is trapped in the jaws of the instrument. Lightly tighten the handle by hand to secure the stem.



3. The extraction rod is screwed into the swivel block.



4. The extraction process is ready to start.



Tel. : +33 (0)2.23.35.53.40 [contact@orthoriginal.bzh](mailto:contact@orthoriginal.bzh)  
Fax : +33 (0)2.23.35.16.35 [www.orthoriginal.bzh](http://www.orthoriginal.bzh)



Orthoriginal, 44C Rue de Bray  
35510 Cesson-Sévigné, France

**Designation**

**Catalogue number**

**Femoral stem extractor**

**0081-01-00**