



ORTH ORIGINAL

Screw extraction set

A simple and efficient solution for extraction of most Torx, hexagonal and cruciform screws




Handle with
AO fitting


TORX tip
T4, T5, T6, T7, T8, T10





ALLEN tip (mm)
1,5-2,0-2,5-3,5



**2 tips for broken screws
+ 1 tip for damaged screw
marks
(for small screws)**



CROSS HEAD tip
Size 1 and Size 2



ORTHORIGINAL
44C rue de Bray
35510 Cesson-Sévigné, FRANCE

Tél : +33 (0)2 23 35 53 40
Fax : +33 (0)2 23 35 16 35
contact@orthoriginal.bzh

